

minirò H-net evo

Fully validatable rotary medical sealer, with sealing parameters control and integrated printer.

The Evolution of one of the most successful Gandus Saldatrici medical sealer.

The basic solution for all sealing requirements of both CSSD and industrial, looking for a medical sealer that fully complies with ISO standards.

User friendly thanks to extended alphanumeric keyboard and bright 4 lines display.

Print data can be stored on the machine with a code and recalled directly from memory.

Integrated printer with 1 line of print.

- **Fully validatable to the EN ISO 11607-2 standards**

Accessories: Front support tray, roller support tray, roll holder complete with cutting device, labels printer, Barcode scanner.

MODELLO	MINIRÒ H-NET EVO
Kind of sealing	Rotary
Sealing width	15 (multi line)
Sealable materials	Pouches/Reels EN 868-5, Tyvek® Pouches, Paper pouches EN 868-4, Gusseted pouches, Laminated AL pouches, PP pouches
ISO 11607-2 Validation	Yes
DIN 58953-7:2008 Conformity	Yes
Printer	Built-in
Interfaces	2 serial ports RS232
Dimensions (mm)	554 x 270 x 222
Weight (kg)	18,5

