

minirò H-seal VP evo

Fully validatable rotary sealer with complete control of sealing parameters.

The basic solution for all sealing requirements of CSSD and industries looking for a medical sealer that fully complies with ISO standards.

Operator's friendly thanks to wide keyboard and bright 2 lines with 26 characters display.

All functioning data can be stored onboard, immediately recalled from memory and easily exportable through USB.

It can be linked with label printer.

- **Fully validatable to the EN ISO 11607-2 standards**

Accessories: Front support tray, roller support tray, roll holder complete with cutting device, labels printer.

MODEL	MINIRO' H-SEAL VP EVO
Kind of sealing	Rotary
Sealing width (mm)	15 (multiline)
Sealing length (mm)	Unlimited
Sealable materials	Pouches/Reels EN 868-5, Tyvek® Pouches, Paper pouches EN 868-4, Gusseted pouches, Laminated AL pouches, PP pouches
ISO 11607-2 Validation	Yes
DIN 58953-7:2008 Directive	Yes
Printer	External
Interfaces	1 RS232 port, 1 USB port
Dimensions (mm)	554 x 270 x 222
Weight (kg)	16

