

minirò poly VP

**Rotary sealer for PE, PP, PP/PA, Tyvek® and heat sealable laminated materials pouches.
 Fully validatable according to ISO 11607-2 rules.**

Traceability system.

External print system available.

Speed 6 mt/min.

Temperature adjustment 10-220°C and machine stop in case of overheating.

Automatic sealing pressure adjustment.

LCD 2 lines, 16 characters.

Materials: PE, PP, laminated, carton head bags.

Accessories: motorized conveyor, PS 654 front support tray, RL 80 roller support tray.

MODELLO	minirò poly VP
Kind of sealing	Rotary – knurled (flat seal option)
Sealing height	9 mm
Free edge	adjustable 0 - 32 mm
Temperature adjustment	10-220°C
Speed	6 m/min (optional 4m/min)
Sealable materials	PE-PP-laminated-carton head bags
ISO 11607-2 Validation	YES
DIN 58953-7:2008 Conformity	YES
Printer	via RS 232 to external label printer
Interface with external connections	1 serial port RS232 - 1 USB port
Dimensions (mm)	590 x 270 x 238
Weight (kg)	18

